

# Small Adult Disposable Manual Resuscitators with Built-In Pressure Monitor

*Single Patient Use Non-Sterile*

*Re-Order AF5000 & AF5100 Series*

**\*\*Non-Latex\*\***

## Small Adult AirFlow Resuscitator

AF 5040 (without Manometer)

AF 5140 (with Manometer)



Bag Volume	1500 ml		
Stroke Volume	780 ML (ASTM Standard Hand)		
Deadspace	10 ml		
Exhalation Resistance	< 4 cm H <sub>2</sub> O @ 50 lpm flow		
Pressure Relief (optional)	40 cm H <sub>2</sub> O		
Attainable Delivery Pressure	≥ 100 cm H <sub>2</sub> O		
Ventilatory Frequency	75 - 80 bpm (NO LOAD)		
Operating Temperature	-18°C to 51°C (0o TO 123°F)		
Storage Environment Limits	-40°C to 60°C (-40o TO 140°F)		
Patient Port Inlet	15mm/22mm		
Average Device Mass without Mask	0.38 kg (0.84 lb)		
Average Length (Basic Model)	33 cm (13")		
Body Mass Use Range	≥ 40 kg (88 lbs)		
TIDAL VOLUME	550 ml	550 ml	550 ml
Rate	12 bpm	15 bpm	20 bpm
O <sub>2</sub> flowrate	15 lpm	15 lpm	15 lpm
F <sub>D</sub> O <sub>2</sub>	99%	99%	99%