

Small Child Disposable Manual Resuscitators with Built-In Pressure Monitor ****Non Latex****

With 40cmH₂O Pressure Limiting System & Pressure Port Adapter
Single Patient Use Non-Sterile ***Re-Order AF4000 & AF4100 Series***

Small Child AirFlow Resuscitator

AF 4040 (without Manometer)

AF 4140 (with Manometer)



Bag Volume	550 ml		
Stroke Volume	> 350 ML (ASTM Standard Hand)		
Deadspace	10 ml with mask or low dead space adapter		
Exhalation Resistance	< 4 cm H ₂ O @ 5 lpm flow		
Pressure Relief (optional)	40 cm H ₂ O		
Attainable Delivery Pressure	40 cm H ₂ O and override		
Ventilatory Frequency	> 80 bpm (NO LOAD)		
Operating Temperature	-18°C to 51°C (0o TO 123°F)		
Storage Environment Limits	-40°C to 60°C (-40o TO 140°F)		
Patient Port Inlet	15mm/22mm		
Average Device Mass without Mask	0.22 kg (0.50 lb)		
Average Length (Basic Model)	25 cm (10")		
Body Mass Use Range	< 22 kg (48 lbs)		
TIDAL VOLUME	300 ml	300 ml	300 ml
Rate	12 bpm	20 bpm	30 bpm
O ₂ flowrate	10 lpm	10 lpm	10 lpm
F _D O ₂	99%	99%	99%