

Adult Disposable Manual Resuscitators with Built-In Pressure Monitor

Single Patient Use Non-Sterile Re-Order AF1000 & AF1100 Series
****Non-Latex****

Adult AirFlow Resuscitator
AF 1040 (without Manometer)
AF 1140 (with Manometer)



Bag Volume	1900 ml		
Stroke Volume	≥ 960 ML (ASTM Standard Hand)		
Deadspace	10 ml		
Exhalation Resistance	< 4 cm H ₂ O @ 50 lpm flow		
Pressure Relief (optional)	40 cm H ₂ O		
Attainable Delivery Pressure	≥ 100 cm H ₂ O		
Ventilatory Frequency	80 - 85 bpm (NO LOAD)		
Operating Temperature	-18°C to 51°C (0o TO 123°F)		
Storage Environment Limits	-40°C to 60°C (-40o TO 140°F)		
Patient Port Inlet	15mm/22mm		
Average Device Mass without Mask	0.46 kg (1.02 lb)		
Average Length (Basic Model)	35 cm (14")		
Body Mass Use Range	≥ 40 kg (88 lbs)		
TIDAL VOLUME	600 ml	600 ml	600 ml
Rate	12 bpm	15 bpm	20 bpm
O ₂ flowrate	15 lpm	15 lpm	15 lpm
F _D O ₂	99%	99%	99%