

## The Ventlab Hyperinflation System

Our Hyperinflation systems are available in 5 sizes with either no manometer or a built-in manometer. You may also choose to include a face mask to complete the system.

The Hyperinflation System with Manometer has a built-in color coded manometer with a 40 cm Pressure Relief Pop-Off valve, Swivel Patient Connection Elbow and a Latex-Free triple pleated bag (1/4 liter bag has two pleats). As with all Ventlab product, the Hyperinflation system is completely *Latex-Free*.



**HS1000 Series** - Hyperinflation System includes standard patient connection elbow, oxygen connection port, manometer port, latex free breathing bag, O<sub>2</sub> tubing and pressure tubing.



**HS 4000 Series** - Hyperinflation System includes swivel patient connection elbow with manometer, latex free breathing bag and O<sub>2</sub> tubing.



Manometer with pop-off valve – closed



Manometer with pop-off valve – open

### Ordering Information

Choose your product by using the table below.

Replace the first red number with the desired size face mask (0 for no mask) and the second for bag size as shown in the example.

**Case Qty: 20 each**

Hyperinflation System - With Manometer	Hyperinflation System - No Manometer
Example: HS4 <sup>2</sup> 5 <sup>1</sup> = Infant Mask and ½ liter Bag	Example: HS1 <sup>2</sup> 5 <sup>0</sup> = Infant Mask and ½ liter Bag
<div style="display: flex; justify-content: space-around;"> <div style="text-align: left;"> <p><b>Masks</b></p> <p>0 = No Mask</p> <p>1 = SM Infant</p> <p>2 = Infant</p> <p>3 = Sm. Child</p> <p>4 = Child</p> <p>5 = Adult (reg)</p> <p>6 = Adult (large)</p> <p>7 = Adult (reg. wide)</p> </div> <div style="text-align: left;"> <p><b>Bags</b></p> <p>4 = ¼ liter</p> <p>5 = ½ liter</p> <p>1 = 1 liter</p> <p>2 = 2 liter</p> <p>3 = 3 liter</p> </div> </div>	<div style="display: flex; justify-content: space-around;"> <div style="text-align: left;"> <p><b>Masks</b></p> <p>0 = No Mask</p> <p>1 = SM Infant</p> <p>2 = Infant</p> <p>3 = Sm. Child</p> <p>4 = Child</p> <p>5 = Adult (reg)</p> <p>6 = Adult (large)</p> <p>7 = Adult (reg. wide)</p> </div> <div style="text-align: left;"> <p><b>Bags</b></p> <p>4 = ¼ liter</p> <p>5 = ½ liter</p> <p>1 = 1 liter</p> <p>2 = 2 liter</p> <p>3 = 3 liter</p> </div> </div>