

VENTLAB INTRODUCES NEW CO₂ DETECTION DEVICE

StatCheck II™

Ventlab is happy and proud to announce a “new and improved” version of our StatCheck™ CO₂ indicator. The new **StatCheck II™** follows the popular convention of “PURPLE to BEIGE to YELLOW” color change. Our previous version changed from White to Purple. Although we offered many clinical and financial advantages over competitive CO₂ detectors, the fact that the majority of devices turn from PURPLE to YELLOW limited the acceptance of StatCheck™ in the marketplace.

Effective immediately, Ventlab’s **StatCheck II™** products will turn from **PURPLE to BEIGE to YELLOW**, instead of White to Purple.



The new and improved **StatCheck II™** not only conforms to the more popular color scheme, but also offers other clinical advantages:

- ✓ Can be pre-attached and ready for instant use (thus improving compliance). It may be left attached to continuously monitor any risk of tube migration during resuscitation.
- ✓ Quick and simple to deploy. To activate, simply turn cap clockwise to tighten, and then remove it.
- ✓ May be used on patients as light as 2.5kg. No need to inventory multiple sizes.
- ✓ Very competitively priced. Reduces the cost of CO₂ detection.
- ✓ Side stream design eliminates chance of vomitus contamination.

CALL VENTLAB @ 800-593-4654 FOR YOUR FREE SAMPLE TODAY!