



NASAL CPAP CANNULA SYSTEM

VENTLAB introduces our new Nasal Cannula System for CPAP and Bi-Level therapy. This single patient use disposable product enhances patient comfort and compliance. Our less obtrusive design allows for better communication between patients and clinicians.



Wider Field of Vision

Multiple Adjustments for Greater Stability and Comfort

Ports for Gas Delivery and Monitoring

Latex free

Comfortable Featherlite Padded Headgear

Soft Slimline Circuit

Comfort Fit Nasal Inserts



Available in Four Sizes

Purified Adapters

From Saving Lives to Improving the Quality of Lives... You can trust VENTLAB

The **VENTLAB** Nasal Cannula CPAP System consists of headgear, padding, circuit, tubing clips and a single nasal cannula. Two pieces of flexible tubing join the cannula to a "Y" shaped coupling for connection to a positive pressure device. The system includes ports for pressure monitoring, gas sampling and O₂ delivery.

Specifications

- 1.7 cmH₂O @ 40 lpm
- 3.3 cmH₂O @ 60 lpm
- 5.6 cmH₂O @ 80 lpm
- 8.3 cmH₂O @ 100 lpm
- Sound Level Pressure 16 db



Ordering Information

- Small Nasal CPAP Cannula System (10 per case)
- Medium Nasal CPAP Cannula System (10 per case)
- Large Nasal CPAP Cannula System (10 per case)
- X-Large Nasal CPAP Cannula System (10 per case)
- Small Headgear w/Clips and Padding (20 per case)
- Large Headgear w/Clips and Padding (20 per case)

Item Number

- CS4000
- CS4001
- CS4002
- CS4003
- CS100
- CS200