



NASAL CPAP CANNULA SYSTEM

VENTLAB introduces our new Nasal Cannula System for CPAP and Bi-Level therapy. This product offers the benefits needed to ensure patient comfort and satisfaction. Our design enhances patient mobility and provides the comfort people demand.

Wider Field of Vision

Multiple Adjustments for Greater Stability and Comfort

Latex Free

Soft Slimline Circuit

Comfortable Featherlite Padded Headgear

Comfort Fit Nasal Inserts

Easy to Clean



Each Kit Contains Four Sizes



VENTLAB understands your need for comfort and system stability. Our Nasal CPAP Cannula System is packaged complete with all sizes of headgear and nasal inserts. It provides the freedom of motion and field of vision your patients are accustomed to. Ventlab's innovative design encourages compliance which helps promote long term success.

From Saving Lives to Improving the Quality of Lives... You can trust VENTLAB



The **VENTLAB** Nasal CPAP Cannula System contains:

- 4 sizes of nasal inserts
- 2 sizes of adjustable headgear
- Patient circuit with swivel connector
- Tubing clips
- Headgear padding
- Convenient nylon carrying case

Specifications

1.7 cmH₂O @ 40 lpm

3.3 cmH₂O @ 60 lpm

5.6 cmH₂O @ 80 lpm

8.3 cmH₂O @ 100 lpm

Sound Level Pressure 16db

Ordering Information

Ventlab Reusable Nasal CPAP Cannula System (6 per Case)

Large Headgear w/clips and padding

Small Headgear w/clips and padding

Small Cannula

Medium Cannula

Large Cannula

X-Large Cannula

Patient Circuit

Tubing Clips

Headgear Padding



Item Number

CS9000

CS9001

CS9002

CS9003

CS9004

CS9005

CS9006

CS9007

CS9008

CS9009