

Ventlab Adult Resuscitator VN1000 Series
 (All Ventlab Resuscitation Bags are latex free)



PARAMETERS	VENTLAB ADULT RESUSCITATORS SPECS.	ASTM STD.
Bag Volume	1900 ml	
Stroke Volume	One hand 780-980 ml	600 ml
Deadspace	10 ml	<30 ml
Max. Cycles	80-85 BPM (no load)	20 bpm
FDO ₂	95% - 98% attainable at 15 lpm O ₂ , 12 bpm, 600 – 800 ml stroke volume (VN 1000MX)	>85% at 600 ml
PEEP Capability	Use any make 19 mm ID PEEP valve	
Pressure Relief	Available	
Operating Temp	0°F to 123°F	0°F to 123°F
Drop Test	Safe at height of 1 meter	1 meter height

PRODUCT NO.	DESCRIPTION	PACKING
VN1000MX	Adult Bag, Mask w/Collapsible Accumulator	10 ea./case
VN 1000MB	Adult Bag, Mask w/Inflatable Bag O ₂ Reservoir	10. ea./case
VN 1000MC	Adult Bag, Mask w/Corr. Tubing O ₂ Reservoir	10 ea./case
VN 1000OX	Adult Bag, w/Collapsible Accumulator, No Mask	10 ea./case
VN1000MXP	Adult Bag, PEEP w/Collapsible Accumulator	10 ea./case