

**Ventlab Infant Resuscitator VN3000 Series**  
(All Ventlab Resuscitation Bags are latex free)



PARAMETERS	VENTLAB INFANT RESUSCITATORS SPECS.	ASTM STD.
Bag Volume	300 ml	
Stroke Volume	One hand 10 – 15 ml	
Deadspace	7 ml with mask or low dead space adpt.	<7 ml
Max. Cycles	70-90 BPM (no load)	20 bpm
FDO <sub>2</sub>	95% - 98% attainable at 15 lpm O <sub>2</sub> , 60 bpm, 10 – 150 ml stroke volume	>85% at 600 ml
PEEP Capability	Use any make 19 mm ID PEEP valve	
Pressure Relief	40 cm H <sub>2</sub> O	
Operating Temp	0°F to 123° F	0° F to 123° F
Drop Test	Safe at height of 1 meter	1 meter height

**ORDERING INFORMATION**

PRODUCT NO.	DESCRIPTION	PACKING
VN 3000MB	Infant Bag, Mask w/Inflatable Bag O <sub>2</sub> Reservoir	10. ea./case
VN 3000MC	Infant Bag, Mask w/Corr. Tubing O <sub>2</sub> Reservoir	10 ea./case