

Ventlab Small Adult Resuscitator VN5000 Series

(All Ventlab Resuscitation Bags are latex free)



PARAMETERS

BAG VOLUME

STROKE VOLUME

DEADSPACE

EXHALATION RESISTANCE

PRESSURE RELIEF(optional)

ATTAINABLE DELIVERY PRESSURE

VENTILATORY FREQUENCY

OPERATING TEMPERATURE

STORAGE ENVIRONMENT LIMITS

PATIENT PORT INLET

AVG. DEVICE MASS w/o MASK

AVG. LENGTH (BASIC MODEL)

BODY MASS USE RANGE

SPECIFICATIONS

1600 ml

780 ml (ASTM Standard Hand)

10 ml

<4 cm H₂O @ 50 lpm flow

40 cm H₂O

≥100 cm H₂O

60-65 bpm (NO LOAD)

-18°C TO 51°C (0° TO 123°F)

-40°C TO 60°C (-40°F TO 140°F)

15 mm/22 mm

0.38 kg (0.84lb.)

29.85 cm (11.75")

≥**40 kg (88 lbs.)**

VN 5000 SERIES DELIVERED OXYGEN CONCENTRATION

Tidal Volume	550 ml	550 ml	550 ml
Rate	12 bpm	15 bpm	20 bpm
Oxygen Flowrate	15 lpm	15 lpm	15 lpm
FDO ₂	99%	99%	99%

ORDERING INFORMATION

Product Number	Description	Qty/Case
VN5000MX	Small Adult Bag, Mask w/Collapsible Accumulator & Built-in Pressure Monitor	10
VN5000MB	Small Adult Bag, Mask w/Inflatable Bag, O2 Reservoir & Built-in Pressure Monitor	10
VN5000MC	Small Adult Bag, Mask w/Corrug. Tubing, O2 Reservoir & Built-in Pressure Monitor	10