

Ventlab Small Adult Resuscitator VN5100 Series

(All Ventlab Resuscitation Bags are latex free)



PARAMETERS

SPECIFICATIONS

BAG VOLUME	1600 ml
STROKE VOLUME	780 ml (ASTM Standard Hand)
DEADSPACE	10 ml
EXHALATION RESISTANCE	<4 cm H ₂ O @ 50 lpm flow
PRESSURE RELIEF(optional)	40 cm H ₂ O
ATTAINABLE DELIVERY PRESSURE	≥100 cm H ₂ O
VENTILATORY FREQUENCY	60-65 bpm (NO LOAD)
OPERATING TEMPERATURE	-18°C TO 51°C (0° TO 123°F)
STORAGE ENVIRONMENT LIMITS	-40°C TO 60°C (-40°F TO 140°F)
PATIENT PORT INLET	15 mm/22 mm
AVG. DEVICE MASS w/o MASK	0.38 kg (0.84lb.)
AVG. LENGTH (BASIC MODEL)	29.85 cm (11.75")
BODY MASS USE RANGE	≥ 40 kg (88 lbs.)

VN 5100 SERIES DELIVERED OXYGEN CONCENTRATION

Tidal Volume	550 ml	550 ml	550 ml
Rate	12 bpm	15 bpm	20 bpm
Oxygen Flowrate	15 lpm	15 lpm	15 lpm
FDO ₂	99%	99%	99%

ORDERING INFORMATION

Product Number	Description	Qty/Case
VN5100MX	Small Adult Bag, Mask w/Collapsible Accumulator & Built-in Pressure Monitor	10
VN5100MB	Small Adult Bag, Mask w/Inflatable Bag, O2 Reservoir & Built-in Pressure Monitor	10
VN5100MC	Small Adult Bag, Mask w/Corrug. Tubing, O2 Reservoir & Built-in Pressure Monitor	10