

The Ventlab Universal Circuit for CPAP and Bi-Level

Reorder part # UC1000 Single Patient Use Latex Free Rx Only



The Ventlab Universal Circuit for CPAP and BiLevel is the most universal single limb circuit on the market. The circuit will interface with flow generators and ventilatory support systems manufactured in the United States, Australia, France and Germany.

The circuit will work with systems that may or may not require:

- An **exhalation port test** to set leak values
- An **external pressure monitoring line** for use with external pressure alarms
- A **luer fitting** to connect to an external alarm package



The **dual ported adapter** features a 22mm ID connection for interface to a variety of full face masks, nasal masks and inserts. It functions as a **pressure monitoring adapter** proximal to the patient, and as an **oxygen entrainment port** or **exhalation port** for non-vented masks.

A standard 22mm OD adapter is included in the circuit package to adapt to non-vented interfaces.